



DryCLO is the next generation in *Helicobacter pylori* (*H. pylori*) testing.

- **Rapid results**
- **No refrigeration**
- **No reagents**

The DryCLO – Rapid Urease test

- Can be stored at room temperature;
- Fast results – positive in as little as 1 minute; negative results within 1 hour;
- Accurate – 99% sensitivity and 97% specificity (as shown in major US studies);
- Is a 'dry test' so there's no incubation or reagents required;
- Simple and convenient – CLIA waived;
- It's cost effective; and
- Good stability with a 2-year shelf life.

5 SIMPLE STEPS



1. Peel back label



2. Place biopsy on test



3. Reseal test

4. Wait up to one hour



5. Read results

H. pylori is a bacteria that lives only in the lining of the stomach and is one of the most common chronic infections in humans. The importance of *H. pylori* was recognised in 1982, when two Australian doctors, Professor Barry Marshall and his colleague Professor Robin Warren, discovered that the germ was almost always present in patients with gastritis (inflammation of the stomach) and ulcers. *H. pylori* causes inflammation of the stomach lining and weakens the natural protection against stomach acid, which may then cause an ulcer. If stomach acid is reduced with medications, an ulcer may heal, but it tends to come back once the medication is stopped. However, if *H. pylori* is treated successfully with antibiotics, an ulcer can be permanently cured. This means that the ulcer medication may no longer be needed.

Distributed in the Philippines by:

**Genesis Lifetech
Marketing**

Antipolo City, Rizal, Philippines

For further information and orders

Genesis Lifetech Marketing
Phone / WhatsApp: +63 908 899 7989
Email: inquiry@genesislifetech.com
www.genesislifetech.com